

mdi E-Funnel is an easy to use pre-sterilized, individually packed, disposable filtration device specially designed to help facilitate and speed up microbial analysis of antibiotics and drugs containing bacteriostats in pharmaceutical API and formulation industries.

E-Funnel has a convenient design with a detachable 100 ml funnel and pluggable bottom for ease of sample collection and transfer to laboratory.

It houses a 0.45µm Edge Hydrophobic Cellulose Nitrate membrane, to improve sensitivity and reliability by ensuring complete flushing of the drug sample, which is critical during bio-burden testing of antibiotics and drugs containing bacteriostats.

Applications

- Product bio-burden testing in pharmaceutical API manufacturing
- Raw material bio-burden testing in pharmaceutical formulation manufacturing

Unique Performance Advantages

- Ready to use
- No sterilization required
- No re-usable stainless steel or polysulfone funnels
- Minimum lab space
- Minimizes chances of extraneous contamination

Ordering Information

Type		Size		Pore Size		XX	XX	Sterility		Pack Size	
Type	Code	Size	Code	Pore Size	Code				Code	Pack Size	Code
E-Funnel	FMEN	47mm	09	0.45µm	02			EO Sterile	2	24	12

Example:

FMEN	09	02	XX	XX	2	12
------	----	----	----	----	---	----



Specifications

Pore Size
0.45µm

Diameter
47 mm

Retention Efficiency
LRV > 7 for *Sr. marcescens*

Sterilization
EO

Water Flow Rates
200ml/min at 250mm Hg Vacuum

Validated for

- Microbial recovery with antibiotic drug samples
- Microbial retention
- Sterility